

PowerHickman™ Central Venous Catheter

Durable Power-Injection CVC



Specially indicated for the power injection of contrast media for CT scans, the PowerLine™ Central Venous Catheter is ideal for central venous pressure monitoring and administration of I.V. fluids, blood products, medications, and blood withdrawals.¹

Power Injectable¹

- Maximum 5 ml/sec flow rate at 300 psi pressure limit enables injection of contrast media for contrast enhanced computed tomography (CECT) scans

Versatility, Strength and Durability

- Single, dual, and triple lumen options help increase the flexibility of your treatment
- Indicated for central venous pressure monitoring¹
- Polyurethane construction adds strength and durability for both short-term and long-term use

Clear Labeling

- Easy to identify purple power injectable lumen
- Depth markings every 1 cm along catheter

The Power of Purple™


Robust luer connectors designed for durability during long-term use

Purple power injectable lumen labeled for easy identification

All extension legs labeled "Cuff" to identify as a tunneled catheter

I.D. tag indicates maximum injection rate



Compatible with Statlock™ CV Plus Stabilization Device¹



SURECUFF™ in-growth cuff promotes tissue growth to secure catheter placement

1 cm depth markings help placement accuracy

Single, dual, and triple lumen catheter configurations



Ordering Information

Lumens	French Size	Tray	Priming Volume (mL)	Average Gravity Flow Rate (mL/hr)	Cuff Location (cm)	Length* (cm)	Overall Length (cm)	Inner Diameter (in.)	Outer Diameter (in.)	Product Code
Single	8F	Intermediate	1.5	6780	5	47	61.7	0.060	0.106	6805810
Single	8F	Microintroducer	1.5	6780	2	47	61.7	0.060	0.106	6802815
Single	8F	Microintroducer	1.5	6780	5	47	61.7	0.060	0.106	6805815
Dual	9.5F	Intermediate	1.3 (each)	3505 (each)	5	47	62.3	0.036	0.126	6805910
Dual	9.5F	Microintroducer	1.3 (each)	3505 (each)	2	47	62.3	0.036	0.126	6802915
Dual	9.5F	Microintroducer	1.3 (each)	3505 (each)	5	47	62.3	0.036	0.126	6805915

* Distance from the most proximal measurement marker to catheter tip.

Tray Components

Intermediate	Microintroducer
<ul style="list-style-type: none"> • Powerhickman™ Polyurethane Catheter with SURECUFF™ Tissue Ingrowth Cuff • Airguard™ Valved Introducer • 18 G Introducer Needle • 45 cm Guidewire (0.035 in.) • Tunneler • Statlock™ Securement Device • Syringe • Dressing, Transparent • Injection Cap(s) • Safety Scalpel • Vein Pick 	<ul style="list-style-type: none"> • Powerhickman™ Polyurethane Catheter with SURECUFF™ Tissue Ingrowth Cuff • Airguard™ Valved Introducer • MicroEZ™ Microintroducer • 21 G Introducer Needle • 70 cm Guidewire (0.018 in.) • 45 cm Guidewire (0.035 in.) • Tunneler • Statlock™ Securement Device • Syringe • Dressing, Transparent • Injection Cap(s) • Safety Scalpel • Vein Pick

References

1. PowerHickman™ Central Venous Catheter Instructions For Use.
2. AirGuard™ Valved Introducer Instructions For Use

Please consult product labels and instructions for use for all indications, contraindications, Hazards and warnings precautions.

벡톤디킨슨코리아(주)
서울특별시 강남구 테헤란로 142, 16층

bd.com

