

# PowerLine™ Central Venous Catheter

Small Diameter Power-Injection CVC



Specially indicated for the power injection of contrast media for CT scans, the PowerLine™ Central Venous Catheter is ideal for central venous pressure monitoring and administration of I.V. fluids, blood products, medications, and blood withdrawals.<sup>1</sup>

## Power Injectable<sup>1</sup>

- Maximum 5 ml/sec flow rate at 300 psi pressure limit enables injection of contrast media for contrast enhanced computed tomography (CECT) scans

## Versatility, Strength and Durability

- Single, dual, and triple lumen options help increase the flexibility of your treatment
- Indicated for central venous pressure monitoring<sup>1</sup>
- Polyurethane construction adds strength and durability for both short-term and long-term use

## Clear Labeling

- Easy to identify purple power injectable lumen
- Depth markings every 1 cm along catheter

# The Power of Purple™

Robust luer connectors designed for durability during long-term use


Purple power injectable lumen labeled for easy identification

All extension legs labeled "Cuff" to identify as a tunneled catheter

I.D. tag indicates maximum injection rate





Compatible with Statlock™ CV Plus Stabilization Device<sup>1</sup>



SURECUFF™ in-growth cuff promotes tissue growth to secure catheter placement

1 cm depth markings help placement accuracy

Single, dual, and triple lumen catheter configurations



## Ordering Information

Lumens	French Size	Tray	Priming Volume (MI)	Average Gravity Flow Rate (mL/hr)	Cuff Location (cm)	Length* (cm)	Overall Length (cm)	Product Code
Single	5F	Microintroducer	0.62	1245	5	50	62.4	6700515
Dual	5F	Microintroducer	0.55 / 0.55	577 / 577	5	50	63.3	6720515
Triple	6F	Microintroducer	0.67 / 0.45 / 0.45	924 / 360 / 360	5	50	65.8	6730615

\* Distance from the most proximal measurement marker to catheter tip.

## Tray Components

### Microintroducer

- Reverse Taper Polyurethane Catheter with SURECUFF™ Tissue Ingrowth Cuff
- MICROEZ™ Microintroducer
- Hydro-Glide™ Insertion Stylet
- 21 G Introducer Needle
- 70 cm Guidewire (0.018 in.)
- Tunneler
- Scalpel
- Attachable Suture Wing
- Dressing, Transparent
- Injection Cap(s)
- Syringe
- Statlock™ Securement Device

**Reference** 1. PowerLine™ Central Venous Catheter Instructions For Use.

Please consult product labels and instructions for use for all indications, contraindications, Hazards and warnings and precautions.

바드코리아(주) 06236 서울시 강남구 테헤란로 142, 15층  
T. 02-2188-2900 F. 02-539-3081

[bd.com/ko-kr](http://bd.com/ko-kr)

